- | | | | |
|-----------|--------------------------------------------------------------------------|------------|---------------------------------------------------------------|
| D | Diagnosewaage
Gebrauchsanweisung 2–10 | TR | Diyagnoz terazisi
Kullanma Talimatı 47–55 |
| GB | Diagnostic scale
Instruction for Use 11–19 | RUS | Диагностические весы
Инструкция по применению 56–65 |
| F | Pèse-personne impédancemètre
Mode d'emploi 20–28 | PL | Waga diagnostyczna
Instrukcja obsługi 66–74 |
| E | Báscula para diagnóstico
Instrucciones para el uso 29–37 | | |
| I | Bilancia diagnostica
Istruzioni per l'uso 38–46 | | |

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Dear customer,

Thank you for choosing one of our products. Our brand stands for well-engineered, quality made products; driven by innovation and safety. Our product range encompasses a wealth of products for health and well-being in the heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage, beauty, baby and air sectors.

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

With kind regards,
Your Beurer team

1. Getting to know your instrument..... 11	8 Taking measurements 15
2. Signs and symbols 11	9. Evaluating the results 16
3. Warnings and safety notes 12	10. Other functions 18
4. Information 13	11. Cleaning and maintaining the device 18
5. Unit description 14	12. Disposal 18
6. Initial use 14	13. What if there are problems? 19
7. Entering user data 14	

Contents

Included in delivery

- Diagnostic scale BF 530
- 2 x 3.0V batteries, CR 2032
- These instructions for use

1. Getting to know your instrument

Functions of the device

This digital diagnostic scale is intended for weighing and providing a diagnosis of your personal fitness data.

It is intended for self-testing in the private domain.

Aside from determining body weight, the scale has the following diagnostic functions that can be used by up to 8 persons:

- BMI,
- Calculating body fat percentage,
- Body water content,
- Muscle percentage,
- Bone mass and
- Basal metabolic rate and active metabolic rate.

The scale also includes the following additional features:

- Switch between kilograms “kg”, pounds “lb” and stones “st”,
- Automatic switch-off function,
- Low battery indicator,
- Automatic user recognition,
- LCD display from P1 to P8

2. Signs and symbols

The following symbols appear in these instructions for use:

WARNING Warning instruction indicating a risk of injury or damage to health.

IMPORTANT Safety note indicating possible damage to the unit/accessory.

Note Note on important information.

3. Warnings and safety notes

Please read these instructions for use carefully and keep them for later use, be sure to also make them accessible to other users and observe the notes they contain.

WARNING

- **The scale may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.**
- Do not use during pregnancy.
- Never step onto the outer edge of the scale with one foot as this could cause the scale to tip!
- Swallowing batteries can be extremely dangerous. Keep the batteries and scale out of the reach of small children. Should a battery be swallowed, seek medical assistance immediately.
- Keep packaging material away from children (risk of suffocation).
- Warning, do not step onto the scale with wet feet or if the surface of the scale is damp – danger of slipping!

Notes on handling batteries

- Replace weak batteries before they discharge completely.
- Always replace all batteries at the same time and use batteries of the same type.
- Batteries must not be recharged, taken apart, thrown into an open fire or short circuited.
- Leaking batteries may damage the device. If you do not intend to use the device for longer periods, remove the batteries from the battery compartment.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Batteries can contain toxins that are harmful to health and the environment. Always dispose of batteries in accordance with applicable legal regulations.
- Do not dispose of batteries with normal household waste.
- Do not throw batteries into open fire. Risk of explosion!

General notes

- The device is suitable only for private use, and not for medical or commercial purposes.
- Please note that measuring tolerances are possible for technical reasons as this scale is not calibrated for use in a professional medical context.
- The scale's maximum capacity is 180 kg (396 lb/28 st). The results for weight and bone mass are displayed in 100 g increments (0.2 lb).
- Measurements of body fat, body water and muscle content are displayed in increments of 0.1%.
- Calorie requirement is indicated in steps of 1 kcal.
- When supplied to the customer, the scale is set to weigh and measure in "kg" and "cm".
- Place the scale on an even, hard surface; a hard surface is crucial to accurate measurements.
- Protect the device from knocks, damp, dust, chemicals, marked temperature fluctuations and nearby sources of heat (ovens, heaters).
- Repairs must only be carried out by Beurer Customer Services or authorised suppliers. Before making a claim, please check the batteries first and replace them if necessary.
- Please contact the specified service address to obtain further information, such as the CE Declaration of Conformity.

Storage and maintenance

The accuracy of the measurements and service life of the device depend on its careful handling:

IMPORTANT

- The device should be cleaned from time to time. Do not use any abrasive cleaning products and never submerge the device in water.
- Make sure that no liquids come into contact with the scale. Never submerge the scale in water. Never rinse it in running water.
- Do not place any objects on the scale when not in use.
- Protect the device from knocks, damp, dust, chemicals, marked temperature fluctuations and nearby sources of heat (ovens, heaters).

- Do not press the buttons violently or with pointed objects.
- Do not expose the scale to high temperatures or strong electromagnetic fields (e.g. mobile telephones).

4. Information

The measurement principle

This scale works on the principle of B.I.A. (bioelectrical impedance analysis). This involves the calculation of body content within seconds, using a current that is not noticeable, is completely non-hazardous and does not pose any risk.

When this measurement of the electrical resistance (impedance) is considered alongside constants and/or individual values (age, height, gender, activity level), the body fat percentage and other variables in the body can be calculated. Muscle tissue and water conduct electricity well and therefore have a lower resistance. In contrast, bones and fat tissue have low conductivity as the fat cells and bones are barely able to conduct the current due to very high levels of resistance.

Please note that the values calculated by the diagnostic scale only represent an approximation of the actual medical analytical values for the body. Only a specialist doctor can precisely establish the body fat, body water, muscle percentage and bone structure using medical methods (e.g. computer tomography).

General tips

- As far as possible, weigh yourself at the same time everyday (ideally in the morning), after having been to the toilet, on an empty stomach and without clothing, in order to obtain results which can be compared.
- Important point regarding the measurement: the calculation of body fat may only be made when bare-foot; the soles of the feet may be lightly dampened where appropriate. Results may be unsatisfactory if the soles of the feet are completely dry or feature high levels of horny skin as this results in insufficient conductivity.
- Stand up straight and still whilst weighing yourself.
- Wait a few hours after any physical exertion to which your body is unaccustomed.
- Wait approx. 15 minutes after getting up so that the water stored in the body can disperse.
- You must remember that only the long-term trend is important. Brief deviations in weight within a few days are normally caused by loss of fluids; however, body water plays an important role in terms of well-being.

Limitations


When calculating the body fat and the additional values, deviating and implausible results may occur in the case of:

- Children under 10 years of age,
- Performance athletes and body-builders,
- Pregnant women,
- Persons with a fever, persons undergoing dialysis, persons demonstrating oedema symptoms and persons suffering from osteoporosis,
- Persons taking cardiovascular medicine (affecting the heart and vascular system),
- Persons taking vasodilative or vasoconstrictive medication,
- Persons with significant anatomical deviations in their legs compared with their overall body height (leg length significantly shorter or longer).

5. Unit description

Overview

1. Electrodes
2. Down button ▼
3. Set button **SET**
4. Up button ▲
5. Display

6. Initial use

Inserting the batteries

If there is one, remove the plastic insulating strip by pulling it from the battery compartment lid and remove the protective film from the batteries, then insert them the correct way round. If the scale does not function, remove the batteries completely and reinsert them.

Changing the weight unit

When supplied to the customer, the scale is set to weigh in “kg”. Switch on the scale by briefly applying pressure to the standing surface with your foot (approx. 1-2 seconds). Wait until “0.0 kg” appears in the display. To the rear of the scale, there is a button where you can select “pound” and “stones” (lb, st).

Positioning the scale

Place the scale on an even, hard surface (no carpet); a hard surface is crucial to accurate measurements.

7. Entering user data

To calculate your body fat percentage and additional body values, you must save the personal user data. The scale has 8 user memory spaces which allow you and members of your family to save and access personal settings, for example.

- **Place the scale on an even, hard surface (no carpet); a hard surface is crucial to accurate measurements.**
- Switch on the scale by briefly applying pressure to the standing surface with your foot (approx. 1-2 seconds). Wait until “0.0 kg” appears in the display.
- Start the setting by pressing the **SET** button. The first memory space flashes in the display.
- Press the ▲/▼ button to select the desired memory space and confirm by pressing the **SET** button.
- You can now make the following settings:

User data	Set values
Users/persons	P-1 to P-8
Height	100 to 220 cm (3' 3.5" to 7' 2.5")
Age	10 to 100 years
Gender	Male (♂), female (♀)
Activity level	1 to 5

Activity levels

The medium-term and long-term perspective is decisive when selecting the activity level.

Activity level	Physical activity
1	None.
2	Low: little and light physical exertion (e.g. going for walks, easy gardening, simple gymnastics).
3	Medium: physical exertion, at minimum 2 to 4 times a week for 30 minutes each.
4	High: physical exertion, at minimum 4 to 6 times a week for 30 minutes each.
5	Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.

- Changing the values: Press the ▲ button or the ▼ button, or press and hold to quickly scroll through the options.
- Confirming the entries: Press the **SET** button.
- Once the values have been set, “0.0 kg” is displayed for 20 seconds.
- Now stand on the scale **with bare feet** whilst “0.0 kg” is displayed.
- Following the weight measurement the BMI, body fat (BF), water, muscle percentage, bone mass, BMR and AMR are displayed

i Note: Your measured weight has now been assigned to your selected user and saved. This is necessary for automatic user detection.

8 Taking measurements

Measuring your weight and carrying out diagnostics

Step onto the scale **with bare feet** and ensure that you are standing still with equal weight distribution and with both legs on the electrodes

i Note: There must be no skin contact between your feet, legs, calves or thighs. Otherwise, the measurement cannot be taken correctly.
The measurement result will be incorrect if the measurement is carried out with socks on.


The scale begins to take measurements straight away. First, the weight is displayed.

When the user stands barefoot on the electrodes, the B.I.A. diagnostics is started automatically.

If a user has been assigned, the user is briefly displayed (e.g. P-); afterwards, the BMI, body fat (BF), water, muscle percentage, bone mass, BMR and AMR are displayed.

The following values appear:

1. Weight in kg
2. BMI
3. Body fat in %
4. Body water in % ≈
5. Muscle percentage in % ↔
6. Bone mass in kg ↔
7. Basal metabolic rate in kcal (BMR)
8. Active metabolic rate in kcal (AMR)

If the user is not recognised, only the weight can be shown and “P--” appears.

Only measuring weight

Now step on the scale wearing **shoes**. Stand still with your weight distributed equally on both feet. The scale begins to take measurements straight away.

The weight is displayed.

i Note: The diagnostic values cannot be calculated as B.I.A. diagnostics is not possible when wearing shoes.

Switching off the scale

The scale switches itself off automatically.

9. Evaluating the results

Body mass index (BMI)

The body mass index (BMI) is a number that is often called upon to evaluate body weight. The figure is calculated from body weight and height. The formula is: $\text{body mass index} = \text{body weight} : \text{height}^2$. The measurement unit for BMI is $[\text{kg}/\text{m}^2]$. According to the BMI, weight is classified for adults (20 years and over) using the following values:

Category		BMI
Underweight	Severely underweight	< 16
	Underweight	16-16.9
	Slightly underweight	17-18.4
Normal weight		18.5-25
Overweight	Overweight	25.1-29.9
Obese (overweight)	Class I obesity	30-34.9
	Class II obesity	35-39.9
	Class III obesity	≥ 40

Body fat percentage

The subsequent body fat percentages serve as a guideline (for further information ask your doctor).

Male

Age	Low	Normal	High	Very high
10-14	<11%	11-16%	16.1-21%	>21.1%
15-19	<12%	12-17%	17.1-22%	>22.1%
20-29	<13%	13-18%	18.1-23%	>23.1%
30-39	<14%	14-19%	19.1-24%	>24.1%
40-49	<15%	15-20%	20.1-25%	>25.1%
50-59	<16%	16-21%	21.1-26%	>26.1%
60-69	<17%	17-22%	22.1-27%	>27.1%
70-100	<18%	18-23%	23.1-28%	>28.1%

Female

Age	Low	Normal	High	Very high
10-14	<16%	16-21%	21.1-26%	>26.1%
15-19	<17%	17-22%	22.1-27%	>27.1%
20-29	<18%	18-23%	23.1-28%	>28.1%
30-39	<19%	19-24%	24.1-29%	>29.1%
40-49	<20%	20-25%	25.1-30%	>30.1%
50-59	<21%	21-26%	26.1-31%	>31.1%
60-69	<22%	22-27%	27.1-32%	>32.1%
70-100	<23%	23-28%	28.1-33%	>33.1%

A lower value will often be determined in the case of athletes. Depending on the sporting activities performed, the training intensity and the physical constitution, the results may still be below the stated standard values. Please note, however, that there may be a risk to health in the case of extremely low values.

Body water content

The body water content in % is normally within the following ranges:

Male

Age	Poor	Good	Very good
10-100	<50%	50-65%	>65%

Female

Age	Poor	Good	Very good
10-100	<45%	45-60%	>60%

Body fat contains relatively little water, therefore body water content may be below the standard values in persons with a high body fat percentage. On the other hand, endurance athletes may exceed the standard values due to having low levels of fat and a high muscle percentage.

The body water calculation performed using this scale is not suitable for drawing medical conclusions, for example in the case of age-related water retention. Consult your doctor where necessary. A high body water content is generally desirable.

Muscle percentage

The muscle percentage is normally within the following ranges:

Male

Age	Low	Normal	High
10-14	<44%	44-57%	>57%
15-19	<43%	43-56%	>56%
20-29	<42%	42-54%	>54%
30-39	<41%	41-52%	>52%
40-49	<40%	40-50%	>50%
50-59	<39%	39-48%	>48%
60-69	<38%	38-47%	>47%
70-100	<37%	37-46%	>46%

Female

Age	Low	Normal	High
10-14	<36%	36-43%	>43%
15-19	<35%	35-41%	>41%
20-29	<34%	34-39%	>39%
30-39	<33%	33-38%	>38%
40-49	<31%	31-36%	>36%
50-59	<29%	29-34%	>34%
60-69	<28%	28-33%	>33%
70-100	<27%	27-32%	>32%

Bone mass

Like the rest of our body, our bones are subject to natural growth, shrinking and ageing processes. The bone mass increases rapidly during childhood and reaches its peak between the ages of 30 and 40. As we age, our bone mass then begins to reduce. You can combat this reduction to an extent with the help of a healthy diet (particularly calcium and vitamin D) and regular physical exercise. The stability of the skeleton can be further increased by means of targeted muscle build-up. Please note that this scale does not identify the calcium content of the bones; instead, it calculates the weight of all the components that make up the bones (organic matter, inorganic matter and water). It is very difficult to influence the bone mass, although it does fluctuate within the scope of influencing factors (weight, height, age, gender). There are no recognised guidelines or recommendations.

IMPORTANT:

Please do not confuse bone mass with bone density.

The bone density can only be determined by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions on changes to bones and bone hardness (e.g. osteoporosis) using this scale.

BMR

The basal metabolic rate (BMR) is the amount of energy required by the body at complete rest in order to maintain its basic functions (e.g. when lying in bed 24 hours a day). This value is primarily dependent on weight, height and age.

It is displayed on the diagnostic scale in the unit kcal/day and calculated using the scientifically recognised Harris-Benedict equation.

This is the amount of energy that is required by your body under all circumstances and must be re-supplied to the body in the form of food. If your energy intake is below this level in the longer term, this can result in damage to health.

AMR

The active metabolic rate (AMR) is the amount of energy that the active body consumes each day. An individual's energy requirement increases when the level of physical activity increases; the diagnostic scale calculates this by means of the activity level (1–5) that has been entered.

To maintain the current weight, the energy that the body uses must be replaced in the form of food and drink.

If less energy is taken in than is used over a long period of time, the body takes the difference from the fat stores with resulting weight loss. However, if the amount of energy taken in exceeds the calculated active metabolic rate (AMR) for a longer period, the body cannot burn off the excess energy. The excess is stored in the body as fat, leading to weight gain.

Temporal context of the results

Please note that only the long-term trend is important. Brief deviations in weight within a few days are normally caused by loss of fluids.

The interpretation of the results is based on changes in total body weight, percentage of body fat, body water and muscle content, as well as on the length of time over which these changes occur.

Rapid changes within the scope of a few days are to be considered separately from medium-term changes (in the scope of weeks) and long-term changes (months).

It can be said as a basic rule that short-term changes in weight are almost entirely changes in water content, whereas medium-term and long-term changes may also relate to the fat percentage and the muscle percentage.

- If your weight drops in the short term but your body fat percentage increases or stays the same, you have only lost water, for example following a training session, visit to the sauna or a diet aimed only at fast weight loss.
- On the other hand, if your weight increases in the medium term and your body fat percentage drops or stays the same, you may have built up valuable muscle mass.
- If your weight and body fat percentage both fall at the same time, then your diet is working – you are losing fat mass.
- Ideally, you should support your diet with physical activity, fitness or strength training. This enables you to increase your muscle percentage in the medium term.
- Body fat, body water and muscle percentages must not be added (muscle tissue also contains components made of body water).

10. Other functions

User assignment

It is possible to assign up to 8 scale users for a diagnostic measurement (barefoot).

For a new measurement, the scale assigns the measurement to the user whose most recently saved measurement falls within +/- 2 kg of this.

Unknown measurements

If measurements cannot be assigned to any user, the display shows “P--”.

Replacing batteries

Your scale is equipped with a low battery indicator. If you operate the scale with low batteries, “Lo” will appear on the display and the scale will automatically switch off. In this case, the batteries must be replaced (2 x 3.0 V, CR 2032).

11. Cleaning and maintaining the device

The device should be cleaned from time to time.

This should be done using a damp cloth and, if necessary, a small amount of detergent.

12. Disposal

Batteries must not be disposed of with household waste. As a consumer, you are required by law to recycle used batteries.

You can recycle your old batteries at public collection points in your community or wherever batteries of the relevant type are sold.


NOTE:

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,

Cd = Battery contains cadmium,

Hg = Battery contains mercury.

For environmental reasons, do not dispose of the scale with the batteries in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point.

Observe the local regulations for material disposal.

Please dispose of the unit in accordance with EC Directive WEEE (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local authorities responsible for waste disposal.

13. What if there are problems?

If the scale encounters an error during measurement, the following is displayed:

Display	Cause	Solution
P--	Unknown measurement (the measurement is beyond the limit of the user assignment and/or it is not possible to assign to a unique user – see chapter 10 "User assignment").	Repeat user assignment. See chapter 7 "Entering user data".
Only the weight is displayed – no diagnostic values are shown.	The proportion of fat is outside the measurable range (less than 3% or greater than 65%).	Please repeat the measurement barefoot, or with the soles of the feet lightly dampened if necessary.
Err	Maximum weight capacity of 180 kg exceeded.	Loads must be limited to 180 kg.
Err or incorrect weight is displayed.	No flat, stable surface.	Place the scale on a flat, stable surface.
Err or incorrect weight is displayed.	Not standing still.	Stand as still as possible.
Incorrect weight is displayed.	Incorrect zero point of scale.	Wait until the scale switches itself off. Activate scale, wait for "0.0 kg" to appear and repeat measurement.
Lo	The batteries in the scale are empty.	Replace the batteries in the scale.

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